



GE Healthcare

Lab Filtration Business  
Cardiff/ Dassel/ Tonglu  
<http://www.gelifesciences.com/>

## EPM2000 Borosilicate Glass Fibre Filters

To whom it may concern:

This is to certify that GE Healthcare Whatman EPM2000 Borosilicate Glass Fibre Filters meet the requirements of:

IS 5182 (Part 4):1999, "Methods for Measurement of Air Pollution – Suspended Particulate Matter".

As referenced under section 7.1, the filter shall have:

A collection efficiency of at least 99% for particles of 0.3µm

Products are manufactured via a Quality System certified to ISO 9001 and tested in accordance with documented quality procedures and approved because of meeting the required specification.

Recommended products (other sizes are available):

1882-866 EPM 2000 Air Sampling Filters, 8 × 10 in sheet (100 pcs)

A handwritten signature in blue ink, appearing to read 'U. Hasting', written over a horizontal line.

U. Hasting  
QA Lead Dassel  
09. Jan. 2017