

Standardowe gilzy z włókna szklanego

Gilzy typu 603g wykonane są z włókien ze szkła borokrzemowego przy zastosowaniu nieorganicznych substancji wiążących. Mogą pracować w temperaturach do 500°C oraz w przypadku rozpuszczalników nieodpowiednich dla gилz celulozowych.

Opcjonalnie dostępna jest także wersja bez lepiszczy.



Rozmiar [mm]*	Grubość ścianki [mm]	Nr katalogowy	Ilość w opak.
Grade 603g (włókna szklane z lepiszczem nieorganicznym)			
10 × 38	1.0	10371103	25
16 × 50	1.0	10371005	25
19 × 90	1.0	10371007	25
22 × 80	1.5	10371011	25
23.8 × 68	1.5	10371114	25
25 × 100	1.5	10371019	25
28 × 60	1.5	10371025	25
30 × 100	1.5	10371036	25
33 × 94	1.5	10371042	25
33 × 100	1.5	10371043	25
33 × 118	1.5	10371045	25
35 × 150	1.5	10371055	25
44 × 230	1.5	10371075	25
Włókno szklane (bez lepiszczy)			
30 × 80	-	2811-308	25

* Średnica wewn. x wysokość zewn.

labsystemservice.pl

Lab-System-Service
ul. Relaksowa 7
70-892 Szczecin
tel. 91 46 223 23, fax 91 46 217 63
e-mail: biuro@s-und-s.pl

GE, the GE Monogram, Anotop, Mini-UniPrep, SPARTAN, Whatman, and Whatman GD/X are trademarks of General Electric Company.
All other third-party trademarks are the property of their respective owner.
© 2019 Lab-System-Service.



Distributor
GE Healthcare

